



Recombinant Arabidopsis thaliana Phospholipase A1-Ildelta (At2g42690)

Product Code	CSB-BP886702DOA
Abbreviation	At2g42690
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9SJI7
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MATTTTSWEE LLGSKNWDTI LDPLDQSLRE LILRCGDFCQ ATYDAFVNDQ NSKYCGASRY GKSSFFDKVM LENASDYEVV NFLYATARVS LPEGLLLQSQ SRDSWDRESN WFGYIAVTSD ERSKALGRRE IYIALRGTSR NYEWNVLGA RPTSADPLLH GPEQDGSGGV VEGTTFSDSDS EDEEGCKVML GWLTIYTSNH PESKFTKLSL RSQLLAKIKE LLLKYKDEKP SIVLTGHSLG ATEAVLAAVD IAENGSSDDV PVTAVFGCP QVGNKEFRDE VMSHKNLIL HVRNTIDLLT RYPGGLLYV DIGINFVIDT KKSPFLSDSR NPGDWHNLQA MLHVVAGWNG KKGEFKLMVK RSIALVKNKSC EFLKAECLVP GSWWVEKNKG LIKNEDEGEWV LAPVEEPPVP EF
Source	Baculovirus
Target Names	At2g42690
Protein Names	Recommended name: Phospholipase A1-Ildelta EC= 3.1.1.-
Expression Region	1-412
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.