



Recombinant Arabidopsis thaliana F-box protein At1g52490 (At1g52490)

Product Code	CSB-MP886765DOA
Abbreviation	At1g52490
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9SSQ2
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MFSV NKRTQA FSVL FFFFSV PKVIRYLNKQ RSLKEEEKRL EEFQ LLLSKL KDVHLPLDLI VEILKKLPTK SLMRFRCVSK PWSFIISKRR DFVESIMARS LRQPPHKLPV FIFHQCDPGT FFTVSSTFSQ STKPKVSIMP GRNHYNARY QYVRGFICCS SSVYDLVTIY NPTTRQCLPL PKIESMVLSP KRHKHCYFGY DHVMNEYKVL AMVND SQELT QTFHVFTLGR DCPQWRKIRG NIDYELISVS RAGVCIDGTI YYVAVRRKDN ENYGELFMMS FDKSERFYH VRTPETLWSP KCTERGLFNH QGKLGCISSN ENNISMWIME NAEKQEWNSI TFGLL EYPGG DFRTFSGITP AGEIFSMCYP FYNMPLCVHY YNVKQKSYRR VEIESTRLKT RKQPRDIIRV FAIHDFVENT MWL
Source	Mammalian cell
Target Names	At1g52495
Protein Names	Recommended name: F-box protein At1g52490
Expression Region	1-423
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.