



Recombinant Escherichia phage Qbeta Capsid protein

Product Code	CSB-EP886842ELC-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9T0R9
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Enterobacteria phage MX1 (Bacteriophage MX1)
Purity	>85% (SDS-PAGE)
Sequence	MAKLQAITLS GIGKNGDVTL NLNPRGVNPT NGVAALSEAG AVPALEKRVT ISVSQPSRNR KNYKVQVKIQ NPTSC TASGT CDPSVTRSAY ADVTFSFTQY STDEERALVR TELKALLADP MLIDAIDNLN PAY
Source	E.coli
Target Names	N/A
Protein Names	Recommended name: Capsid protein Short name= CPA Alternative name(s): Coat protein
Expression Region	1-133
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.