



Recombinant Sphingomyelin synthase-related 2 (W07E6.3)

Product Code	CSB-MP886966CXY
Abbreviation	W07E6.3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9TYV2
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MCKSCENYGKFFFAMLKICIDNSSNFAMKKDAENENLLHAYEHDEQSSKYFIK PDSRYSI DDDISVKIDGKQLEVRKWPTLLATAMVGVGWSNEVALAWVHERVPDDYHPL PDLFFSHF PEVKERLKIFGILGKFEKSEKKTSEFSSPIRGAIRIAEYIMMILLISALLVMFTHQH RWI VIRRVFFCIAMAYSFRALCVTIFQVPVPSINTYCAPKSNSSLELVAGR VVKMFW SAGIEQ LRPRELCGDLIVSGHTLTIFTAFLVFKTYAPQRLQPLSHIYHVLAFTALFSILLAR KHYM IDIVLGYTVSTRIFMEYHALAASYHNRTFETNPLAWSFWSFFIFIPEC DAPANMH NHLLL YNRSTSSKNVSTLKKSRRSFE
Source	Mammalian cell
Target Names	W07E6.3
Protein Names	Recommended name: Sphingomyelin synthase-related 2
Expression Region	1-381
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.