



Recombinant MO25-like protein 3 (mop-25.3)

Product Code	CSB-EP886967CXY-B
Abbreviation	mop-25.3
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q9TZM2
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MFITWCSND APDSSVPRTR SKTKFTP KVR DLLFHC RYEI VDIVDEGYDM PASIGFYNDI IREFVTDDFC LTSLFEDNGQ DSNNIQHRND GCIWSIFDRL MNTNKFRDFD VIQGTFDTLQ IIFFTNHESA NNFIKNNLPR FMQTLHKLIA CSNFFIQAKS FKFLNELFTA QTNYETRSLW MAEPAFIKLV VLAIQSNKHA VRSRAVSILE IFIRNPRNSP EVHEFIGRNR NVLIAFFFN S APIHYYQGSP NEKEDAQYAR MAYKLLNWDM QRPFTQEQLQ DFEEGWTHQK KMREEQLVRT CFHDNNPRLP KVNHVYRTRI IPNQQFLREP LKFGNPFRQ
Source	E.coli
Target Names	mop-25.3
Protein Names	Recommended name: MO25-like protein 3
Expression Region	1-339
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>