



Recombinant Glycogen synthase kinase-3 (gsk-3)

Product Code	CSB-YP886973CXY
Abbreviation	gsk-3
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q9U2Q9
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	<p>MNKQLLSCSL KSGKQVTMVV ASVATDGV DQ QVEISYYDQK VIGNGSFGVV FLAKLSTTNE MVAIKKVLQD KRFKNRELQI MRKLNHPNIV KLKYFFYSSG EKKDELYLNL ILEYVPETVY RVARHYSKQR QQIPMIYVKL YMYQLLRSLA YIHSIGICHR DIKPQNLLID PESGVLKLCD FGSAYLVRN EPNVSYICSR YYRAPELIFG ATNYTNSIDV WSAGTVMAEL LLGQPIFP GD SGVDQLVEII KVLGTP TREQ IQSMNP NYKE FKFPQIKAHP WNKVFRVH TP AEIDLISKI IEYTPTS RPT PQAACQH AFF DELRNPDARL PSGRPLPTLE MDGPMGTGEV STTSGDVAGP SA</p>
Source	Yeast
Target Names	gsk-3
Protein Names	Recommended name: Glycogen synthase kinase-3 EC= 2.7.11.26
Expression Region	1-362
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>