



Recombinant Human Melanoma-associated antigen C2 (MAGEC2)

Product Code	CSB-MP887003HU
Abbreviation	MAGEC2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UBF1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MPPVPGVPFR NVDNDSPTSV ELEDWVDAQH PTDEEEEEAS SASSTLYLVF SPSSFSTSSS LILGGPEEEE VPSGVIPNLT ESIPSSPPQG PPQGPSQSPL SSCCSSFSWS SFSEESSSQK GEDTGTCQGL PDESSSFTYT LDEKVAELVE FLLLKYEAE PVTEAEMLMV VIKYKDYFPV ILKRAREFME LLFGLALIEV GPDHFCVFAN TVGLTDEGSD DEGMpensLL IILSVIFIK GNCASEEVIW EVLNAVGVYA GREHFVYGE RELLTKVWVQ GHYLEYREVP HSSPPYYEFL WGPRAHSESI KKKVLEFLAK LNNTVPSSFP SWYKDALKDV EERVQATIDT ADDATVMASE SLSVMSSNVS FSE
Source	Mammalian cell
Target Names	MAGEC2
Protein Names	Recommended name: Melanoma-associated antigen C2 Alternative name(s): Cancer/testis antigen 10 Short name= CT10 Hepatocellular carcinoma-associated antigen 587 MAGE-C2 antigen MAGE-E1 antigen
Expression Region	1-373
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene is a member of the MAGEC gene family. It is not expressed in normal tissues, except for testis, and is expressed in tumors of various histological types. This gene and the other MAGEC genes are clustered on chromosome Xq26-q27.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.