



Recombinant Human Protein DBF4 homolog A (DBF4)

Product Code	CSB-BP887010HU
Abbreviation	DBF4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UBU7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MNSGAMRIHS KGHFQGGIQV KNEKNRPSLK SLKTDNRPEK SKCKPLWGVK FYLDLPSVTI SEKLQKDIKD LGGRVEEFLS KDISYLISNK KEAKFAQTLG RISPVPSPES AYTAETTSPH PSHDGSSFKS PDTVCLSRGK LLVEKAIKDH DFIPSNSILS NALSWGKIL HIDDIRYYIE QKKKELYLLK KSSTSVRDGG KRVGSGAQKT RTGRLKKPFV KVEDMSQLYR PFYLQLTNMP FINYSIQKPC SPFDVDKPSS MQKQTQVKLR IQTDGDKYGG TSIQLQLKEK KKKGYCECCL QKYEDLETHL LSEQHRNFAQ SNQYQVDDI VSKLVDFDVE YEKDTPKKKR IKYSVGSLSLSP VSASVLKKTE QKEKVELQHI SQKDCQEDDT TVKEQNFLYK ETQETEKLL FISEPIPHPS NELRGLNEKM SNKCSMLSTA EDDIRQNFTQ LPLHKNKQEC ILDISEHTLS ENDLEELRVD HYKCNIQASV HVSDFSTDNS GSQPKQKSDT VLFPAKDLKE KDLHSIFTHD SGLITINSSQ EHLTVQAKAP FHTPPEEPNE CDFKNMDSLP SGKIHRKVKI ILGRNRKENL EPNAEFDKRT EFITQEENRI CSSPVQSLLD LFQTSEEKSE FLGFTSYTEK SGICNVLDIW EEENSDNLLT AFFSSPSTST FTGF
Source	Baculovirus
Target Names	DBF4
Protein Names	Recommended name: Protein DBF4 homolog A Alternative name(s): Activator of S phase kinase Chiffon homolog A DBF4-type zinc finger-containing protein 1
Expression Region	1-674
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.