



Recombinant Human Lipid transferase CIDEB (CIDEB)

Product Code	CSB-BP887032HU
Abbreviation	CIDEB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UHD4
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MEYLSALNPS DLLRSVSNIS SEFGRRVWTS APPPQRPFVRV CDHKRTIRKG LTAATRQELL AKALETLLL N GVLTLVLEED GTAVDSEDFQ QLLEDDTCLM VLQSGQSWSP TRSGVLSYGL GRERP KHSKD IARFTFDVYK QNPRDLFGSL NVKATFYGLY SMSCDFQGLG PKKVLRELLR WTSTLLQGLG HMLLGISSTL RHAVEGAEQW QKGR LHSY
Source	Baculovirus
Target Names	CIDEB
Protein Names	Recommended name: Cell death activator CIDE-B Alternative name(s): Cell death-inducing DFFA-like effector B
Expression Region	1-219
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.