



Recombinant Human Nuclear fragile X mental retardation-interacting protein 1 (NUFIP1)

Product Code	CSB-BP887039HU
Abbreviation	NUFIP1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UHK0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAEPTSDFET PIGWHASPEL TPTLGPLSDT APPRDSWMFW AMLPPPPPL TSSLPAAGSK PSSESQPPME AQSLPGAPP FDAQILPGAQ PPFDAQSPLD SQQPSPGQPW NFHASTSWYW RQSSDRFPRH QKSFNPAVKN SYYPRKYDAK FTDFSLPPSR KQKKKKRKEP VFHFFCDTCD RGFKNQEKYD KHMSEHTKCP ELDCSFTAHE KIVQFHWRNM HAPGMKKIKL DTPEEIARWR EERRKNYPTL ANIERKKKLE KEKEKRGAVL TTTQYGKMKG MSRHSQMAKI RSPGKNHKWK NDNSRQRAVT GSGSHLCDLK LEGPPEANAD PLGVLINSDS ESDKEEKQHV SVIPKEVTPA LC SLMSSYGS LSGSESEPEE TPIKTEADV AENQVLDSSA PKSPSQDVKA TVRNFSEAKS ENRKKSF EKT NPKRKKDYHN YQTLFEP RTH HPYLLEMLLA PDIRHERNVI LQCVRYIIKK DFFGLDTNSA KSKDV
Source	Baculovirus
Target Names	NUFIP1
Protein Names	Recommended name: Nuclear fragile X mental retardation-interacting protein 1 Alternative name(s): Nuclear FMRP-interacting protein 1
Expression Region	1-495
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.