



# Recombinant Human Nuclear fragile X mental retardation-interacting protein 1 (NUFIP1)

<b>Product Code</b>	CSB-EP887039HU-B
<b>Abbreviation</b>	NUFIP1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9UHK0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAEPTSDFET PIGWHASPEL TPTLGPLSDT APPRDSWMFW AMLPPPPPL TSSLPAAGSK PSSESQPPME AQSLPGAPP FDAQILPGAQ PPFDAQSPLD SQQPSPGQPW NFHASTSWYW RQSSDRFPRH QKSFNPAVKN SYYPRKYDAK FTDFSLPPSR KQKKKKRKEP VFHFFCDTCD RGFKNQEKYD KHMSEHTKCP ELDCSFTAHE KIVQFHWRNM HAPGMKKIKL DTPEEIARWR EERRKNYPTL ANIERKKKLE KEKEKRGAVL TTTQYGKMKG MSRHSQMAKI RSPGKNHKWK NDNSRQRAVT GSGSHLCDLK LEGPPEANAD PLGVLINSDS ESDKEEKQHV SVIPKEVTPA LC SLMSSYGS LSGSESEPEE TPIKTEADV AENQVLDSSA PKSPSQDVKA TVRNFSEAKS ENRKKSF EKT NPKRKKDYHN YQTLFEP RTH HPYLLEMLLA PDIRHERNVI LQCVRYIIKK DFFGLDTNSA KSKDV
<b>Source</b>	E.coli
<b>Target Names</b>	NUFIP1
<b>Protein Names</b>	Recommended name: Nuclear fragile X mental retardation-interacting protein 1 Alternative name(s): Nuclear FMRP-interacting protein 1
<b>Expression Region</b>	1-495
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.