



Recombinant Human Zinc finger HIT domain-containing protein 2 (ZNHIT2)

Product Code	CSB-EP887046HU-B
Abbreviation	ZNHIT2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UHR6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MEPAGPCGFC PAGEVQPARY TCPRCNAPYC SLRCYRTHGT CAENFYRDQV LGELRGCSAP PSRLASALRR LRQQRETEDE PGEAGLSSGP APGGLSGLWE RLAPGEKAAF ERLSRGEAG RLLPPWRPWW WNRGAGPQLL EELDNAPGSD AAELAPAR TPPDSVKDAS AAEPAAAERV LGDVPGACTP VVPTRIPAIV SLSRGPVSPL VRFQLPNVLF AYAHTLALYH GGDDALLSDF CATLLGVSGA LGAQQVFASA EEALQAAAHV LEAGEHPPGP LGTRGAMHEV ARILLGEGPT NQKGYTLAAL GDLAQTLGRA RKQAVAREER DHLYRARKKC QFLLAWTNEN EAALTPLALD CARAHQAHAV VAE EVAALTG ELERLWGGPV PPAPRTLIEE LPS
Source	E.coli
Target Names	ZNHIT2
Protein Names	Recommended name: Zinc finger HIT domain-containing protein 2 Alternative name(s): Protein FON
Expression Region	1-403
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.