



Recombinant Human Ena/VASP-like protein (EVL)

Product Code	CSB-BP887052HU
Abbreviation	EVL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UI08
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSEQSICQAR ASVMVYDDTS KKWVPIKPGQ QGFSRINIYH NTASNTFRVV GVKLQDQQVV INYSIVKGLK YNQATPTFHQ WRDARQVYGL NFASKEEATT FSNAMLFALN IMNSQEGGSPS SQRQVQNGPS PDEMIDIQRRQ VMEQHQQQRQ ESLERRTSAT GPILPPGHPS SAASAPVSCS GPPPPPPPPV PPPPTGATPP PPPPLPAGGA QGSSHDESSM SGLAAAIAGA KLRRVQRPED ASGGSSPSGT SKSDANRASS GGGGGGLMEE MNKLLAKRRK AASQSDKPAE KKEDESQMED PSTSPSPGTR AASQPPNSSE AGRKPWERSN SVEKPVSSIL SRTPSVAKSP EAKSPLQSQP HSRMKPAGSV NDMALDAFDL DRMKQEILEE VVRELHKVKE EIIDAIRQEL SGISTT
Source	Baculovirus
Target Names	EVL
Protein Names	Recommended name: Ena/VASP-like protein Alternative name(s): Ena/vasodilator-stimulated phosphoprotein-like
Expression Region	1-416
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.