



# Recombinant Human La-related protein 6 (LARP6)

<b>Product Code</b>	CSB-BP887112HU
<b>Abbreviation</b>	LARP6
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9BRS8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	AQSGGEARP GPKTAVQIRV AIQEAEDVDE LEDEEEGAET RGAGDPARYL SPGWGSASEE EPSRGHSGTT ASGGENERED LEQEWKPPDE ELIKKLVLDQI EFYFSDENLE KDAFLLKHVR RNKLGYSVK LLTSFKVKH LTRDWRRTAH ALKYSVLEL NEDHRKVRRT TPVPLFPNEN LPSKMLLVYD LYLSPKLWAL ATPQKNGRVQ EKVMEHLLKL FGTFGVISSV RILKPGRELP PDIRRISSRY SQVGTQECAI VEFEEVEAAI KAHEFMITES QGKENMKAVL IGMKPPKKK AKDKNHDEEP TASIHLNKSL NKRVEELQYM GDESSANSSS DPESNPTSPM AGRRHAATNK LSPSGHQQLF LSPNASPCTS PWSSPLAQRK GVSRSPLAE EGRINCSTSP EIFRKCMDYS SDSSVTPSGS PWVRRRRQAE MGTQEKSPGT SPLLSRKMQT ADGLPVGVL RPRGPDNTRG FHGHERSRAC V
<b>Source</b>	Baculovirus
<b>Target Names</b>	LARP6
<b>Protein Names</b>	Recommended name: La-related protein 6 Alternative name(s): Acheron Short name= Achn La ribonucleoprotein domain family member 6
<b>Expression Region</b>	2-491
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.