



Recombinant Human La-related protein 6 (LARP6)

Product Code	CSB-EP887112HU-B
Abbreviation	LARP6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BRS8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	AQSGGEARP GPKTAVQIRV AIQEAEDVDE LEDEEEGAET RGAGDPARYL SPGWGSASEE EPSRGHSGTT ASGGENERED LEQEWKPPDE ELIKKLVLDQI EFYFSDENLE KDAFLLKHVR RNKLGYSVK LLTSFKVKH LTRDWRTTAH ALKYSVVLEL NEDHRKVRRT TPVPLFPNEN LPSKMLLVYD LYLSPKLWAL ATPQKNGRVQ EKVMEHLLKL FGTFGVISSV RILKPGRELP PDIRRISSRY SQVGTQECAI VEFEEVEAAI KAHEFMITES QGKENMKAVL IGMKPPKKK AKDKNHDEEP TASIHLNKSL NKRVEELQYM GDESSANSSS DPESNPTSPM AGRRAATNK LSPSGHQQLF LSPNASPCTS PWSSPLAQRK GVSRSPLAE EGRINCSTSP EIFRKCMDYS SDSSVTPSGS PWVRRRRQAE MGTQEKSPGT SPLLSRKMQT ADGLPVGVLRLPRGPDNTRG FHGHERSRAC V
Source	E.coli
Target Names	LARP6
Protein Names	Recommended name: La-related protein 6 Alternative name(s): Acheron Short name= Achn La ribonucleoprotein domain family member 6
Expression Region	2-491
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.