



# Recombinant Human Ubiquitin-associated domain-containing protein 1 (UBAC1)

<b>Product Code</b>	CSB-YP887116HU
<b>Abbreviation</b>	UBAC1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9BSL1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MFVQEEKIFA GKVLRLHICA SDGAEWLEEA TEDTSVEKLIK ERCLKHCAHG SLEDPKSITH HKLIHAASER VLSDARTILE ENIQDQDVLL LIKKRAPSP PKMADVSAEE KKKQDQKAPD KEAILRATAN LPSYNMDRAA VQTNMRDFQT ELRKILVSLI EVAQKLLALN PDAVELFKKA NAMLDEDEDE RVDEAALRQL TEMGFPENRA TKALQLNHMS VPQAMEWLIE HAEDPTIDTP LPGQAPPEAE GATAAASEAA AGASATDEEA RDELTEIFKK IRRKREFRAD ARAVISLMEM GFDEKEVIDA LRVNNNQNA ACEWLLGDRK PSPEELDKGI DPDSPLFQAI LDNPVVQLGL TNPKTLLAFE DMLNPLNST QWMNDPETGP VMLQISRIFQ TLNRT
<b>Source</b>	Yeast
<b>Target Names</b>	UBAC1
<b>Protein Names</b>	Recommended name: Ubiquitin-associated domain-containing protein 1 Short name= UBA domain-containing protein 1 Alternative name(s): E3 ubiquitin-protein ligase subunit KPC2 Glioblastoma cell differentiation-related protein 1 Kip1 u
<b>Expression Region</b>	1-405
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.