



Recombinant Human THUMP domain-containing protein 3 (THUMPD3)

Product Code	CSB-EP887131HU-B
Abbreviation	THUMPD3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BV44
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MCDIEEATNQ LLDVNLHENQ KSVQVTESDL GSESELLVTI GATVPTGFEQ TAADEVREKL GSSCKISRDR GKIYFVISVE SLAQVHCLRS VDNLFVVVQE FQDYQFKQTK EEVLKDFEDL AGKLPWSNPL KVKINASF KKKAKRKKIN QNSSKEKINN GQEVKIDQRN VKKEFTSHAL DSHILDYYEN PAIKEDVSTL IGDDLASCKD ETDESSKEET EPQVLKFRVT CNRAGEKHCF TSNEAARDFG GAVQDYFKWK ADMTNFDVEV LLNIHDNEVI VGIALTEESL HRRNITHFGP TTLRSTLAYG MLRLCDPLPY DIIVDPMCGT GAIEGATE WSDCFHIAGD NNPLAVNRAA NNIASLLTKS QIKEGKPSWG LPIDAVQWDI CNLPLRTGSV DIIVTDLPFG KRMGSKKRNW NLYPACLREM SRVCTPTTGR AVLLTQDTKC FTKALSGMRH VWRKVDTVWV NVGGLRAAVY VLIRTPQAFV HPSEQDGERG TLWQCKE
Source	E.coli
Target Names	THUMPD3
Protein Names	Recommended name: THUMP domain-containing protein 3
Expression Region	1-507
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.