



Recombinant Human Sister chromatid cohesion protein DCC1 (DSCC1)

Product Code	CSB-BP887133HU
Abbreviation	DSCC1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BVC3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MKRTRDEVDA TLQIAKLNAE ELLPAVHCLG FGPGASGAAA GDFCLLELEP TLCQQLEDGH SLVIRGDKDE QAVLCSKDKT YDLKIADTSN MLLFIPGCKT PDQLKKEDSH CNIHTEIFG FSNYWELRR RRPKLLKLLK LLMENPYEGP DSQKEKDSNS SKYTTEDLLD QIQASEEEIM TQLQVLNACK IGGYWRILEF DYEMKLLNHV TQLVDSESWF FGKVPLNTCL QELGPLEPEE MIEHCLKCYG KKYVDEGEVY FELDADKICR AAARMLLQNA VKFNLAEFQE VWQQSVPEGM VTSLDQLKGL ALVDRHSRPE IIFLLKVDDL PEDNQERFNS LFSLREKWTE EDIAPYIQDL CGEKQTIGAL LTKYSHSSMQ NGVKVYNSRR PIS
Source	Baculovirus
Target Names	DSCC1
Protein Names	Recommended name: Sister chromatid cohesion protein DCC1 Alternative name(s): Defective in sister chromatid cohesion protein 1 homolog
Expression Region	1-393
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.