



Recombinant Human Complement C1q tumor necrosis factor-related protein 4 (C1QTNF4)

Product Code	CSB-BP887148HU
Abbreviation	C1QTNF4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BXJ3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	LGPT PPGGSSELRS AFSAARTTPL EGTSEMAVTF DKVYVNIGGD FDVATGQFRC RVPGAYFFSF TAGKAPHKSL SVMLVRNRDE VQALAFDEQR RPGARRAASQ SAMLQLDYGD TVWLR LHGAP QYALGAPGAT FSGYLVYADA DADAPARGPP APPEPRSAFS AARTRSLVGS DAGPGPRHQP LAFDTEFVNI GGDFDAAAGV FRCRLPGAYF FSFTLGKLPK KTL SVKLMKN RDEVQAMIYD DGASRRREM QSVM LALRR GDA VWLLSHD HDGYGAYS NH GK YITFSGFL VY PDLAPAAP PGLGASELL
Source	Baculovirus
Target Names	C1QTNF4
Protein Names	Recommended name: Complement C1q tumor necrosis factor-related protein 4
Expression Region	17-329
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.