



Recombinant Human FSD1-like protein (FSD1L)

Product Code	CSB-MP887149HU
Abbreviation	FSD1L
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BXM9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MDSQKYCFKE NENVTVDKAC FLISNITIGP ESINLQQEAL QRIISTLANK NDEIQNFIDT LHHTLKG VQE NSSNILSELD EEFDLSYLIL DEVKESMINC IKQEQARKSQ ELQSQISQCN NALENSEELL EFATRSLDIK EPEEFSKAAR QIKDRVTMAS AFRLSLKPKV SDNMTHLMVD FSQERQMLQT LKFLPVPKAP EIDPVECLVA DNSVTVAWRM PEEDNKIDHF ILEHRKTNFD GLPRVKDERC WEIIDNIKGT EYTL SGLKFD SKYMNFRVRA CNKAVAGEYS DPVTLET KAL NFNLDNSSSH LNLKVEDTCV EWDPTGGKGQ ESKIKGKENK GRSGTPSPKR TSVGSRP PAV RGSRRDFTGE SYTVLGDTAI ESGQHYWEVK AQKDCSYSV GVAYKTLGKF DQLGKTNTSW CIHVNNWLQN TFAAKHNNKV KALDVTVP EK IGVFCDFDGG QLSFYDANSK QLLYSFKTKF TQPVLP GFMV WCGGLSLSTG MQVPSAVRTL QKSENGMTGS ASSLN NVVTQ
Source	Mammalian cell
Target Names	FSD1L
Protein Names	Recommended name: FSD1-like protein Alternative name(s): Coiled-coil domain-containing protein 10 FSD1 N-terminal-like protein
Expression Region	1-530
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.