



Recombinant Human Testis-expressed sequence 13B protein (TEX13B)

Product Code	CSB-MP887152HU
Abbreviation	TEX13B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BXU2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MALRPEDPSS GFRHGNVAVF IIEKMARHTK GPEFYFENIS LSWEEVEDKL RAILEDSEVP SEVKEACTWG SLALGVRFAH RQGQLQNRV QWLQGFALH RSAALVLASN LTELKEQQEM ECNEATFQLQ LTETSLAEVQ RERDMLRWKL FHAELAPPQG QGQATVFPGL ATAGGDWTEG AGEQEKEAVA AAGAAGGKGE ERYAEAGPAP AEVLQQLGGG FRQPLGAIVA GKLHLCGAEG ERSQVSTNSH VCLLWAWVHS LTGASSCPAP YLIHILIPMP FVRLLSHTQY TPFTSKGHRT GSNSDAFQLG GL
Source	Mammalian cell
Target Names	TEX13B
Protein Names	Recommended name: Testis-expressed sequence 13B protein
Expression Region	1-312
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene is similar to a mouse gene that is expressed in the testis.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.