



# Recombinant Human Testis-expressed sequence 13B protein (TEX13B)

<b>Product Code</b>	CSB-YP887152HU
<b>Abbreviation</b>	TEX13B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9BXU2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MALRPEDPSS GFRHGNVAV IIEKMARHTK GPEFYFENIS LSWEEVEDKL RAILEDSEVP SEVKEACTWG SLALGVRFAH RQGQLQNRV QWLQGFALH RSAALVLASN LTELKEQQEM ECNEATFQLQ LTETSLAEVQ RERDMLRWKL FHAELAPPQG QGQATVFPGL ATAGGDWTEG AGEQEKEAVA AAGAAGGKGE ERYAEAGPAP AEVLQQLGGG FRQPLGAIVA GKLHLCGAEG ERSQVSTNSH VCLLWAWVHS LTGASSCPAP YLIHILIPMP FVRLLSHTQY TPFTSKGHRT GSNSDAFQLG GL
<b>Source</b>	Yeast
<b>Target Names</b>	TEX13B
<b>Protein Names</b>	Recommended name: Testis-expressed sequence 13B protein
<b>Expression Region</b>	1-312
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This gene is similar to a mouse gene that is expressed in the testis.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.