



# Recombinant Human R-spondin-3 (RSPO3)

<b>Product Code</b>	CSB-BP887154HU
<b>Abbreviation</b>	RSPO3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9BXY4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	QNASRGRRQ RRMHPNVSQG CQGGCATCSD YNGCLSCKPR LFFALERIGM KQIGVCLSSC PSGYYGTRYP DINKCTKCKA DCDTCFNKNF CTKCKSGFYL HLGKCLDNCP EGGLEANNHTM ECVSIVHCEV SEWNPWSPCT KKGKTCGFKR GTETRVREII QHPSAKGNLC PPTNETRKCT VQRKKCQKGE RGKKGREKR KKPNKGESKE AIPDSKSLES SKEIPEQREN KQQQKKRKVQ DKQKSVSVST VH
<b>Source</b>	Baculovirus
<b>Target Names</b>	RSPO3
<b>Protein Names</b>	Recommended name: R-spondin-3 Alternative name(s): Protein with TSP type-1 repeat Short name= hPWTSR Roof plate-specific spondin-3 Short name= hRspo3 Thrombospondin type-1 domain-containing protein 2
<b>Expression Region</b>	22-272
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This gene encodes a member of the thrombospondin type 1 repeat supergene family. In addition, the protein contains a furin-like cysteine-rich region. Furin-like repeat domains have been found in a variety of eukaryotic proteins involved in the mechanism of signal transduction by receptor tyrosine kinases.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.