



Recombinant Human UPF0549 protein C20orf43 (C20orf43)

Product Code	CSB-MP887156HU
Abbreviation	C20orf43
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BY42
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MGCDGGTIPK RHELVKGPKK VEKVDKDAEL VAQWNYCTLS QEILRRPIVA CELGRLYNKD AVIEFLLDKS AEKALGKAAS HIKSIKNVTE LKLSDNPAWE GDKGNTKGDK HDDLQRARFI CPVVGLEMNG RHRFCFLRCC GCVFSERALK EIKAEVCHTC GAAFQEDDVI MLNGTKEDVD VLKTRMEERR LRAKLEKTK KPAAESVSK PDVSEEAPGP SKVKTGKPEE ASLDSREKKT NLAPKSTAMN ESSGKAGKP PCGATKRSIA DSEESEAYKS LFTTHSSAKR SKEESAHWVT HTSYCF
Source	Mammalian cell
Target Names	RTFDC1
Protein Names	Recommended name: UPF0549 protein C20orf43
Expression Region	1-306
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.