



Recombinant Human Fibroblast growth factor-binding protein 2 (FGFBP2)

Product Code	CSB-BP887162HU
Abbreviation	FGFBP2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BYJ0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	Q APRQKQGSTG EEFHFQTGGR DSCTMRPSSL GQGAGEVWLR VDCRNTDQTY WCEYRGQPSM CQAFADPKP YWNQALQELR RLHHACQGAP VLRPSVCREA GPQAHMQQVT SSLKGSPEPN QQPEAGTPSL RPKATVKLTE ATQLGKDSME ELGKAKPTTR PTAKPTQPGP RPGGNEEAKK KAWEHWCWKPFL QALCAFLISF FRG
Source	Baculovirus
Target Names	FGFBP2
Protein Names	Recommended name: Fibroblast growth factor-binding protein 2 Short name= FGF-BP2 Short name= FGF-binding protein 2 Short name= FGFBP-2 Alternative name(s): 37 kDa killer-specific secretory protein Short name= Ksp37
Expression Region	20-223
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.