



# Recombinant Human Tripartite motif-containing protein 15 (TRIM15)

<b>Product Code</b>	CSB-MP887181HU
<b>Abbreviation</b>	TRIM15
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9C019
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MPATPSLKVV HELPACTLCA GPLEDAVTIP CGHTFCRLCL PALSQMGAQS SGKILLCPLC QEEEAETPM APVPLGPLGE TYCEEHGEKI YFFCENDAEF LCVFCREGPT HQAHTVGFLD EAIQPYRDRL RSRLEALSTE RDEIEDVKCQ EDQKLQVLLT QIESKKHQVE TAFERLQQEL EQQRCLLLAR LRELEQQIWK ERDEYITKVS EEVTRLGAQV KELEEKCQP ASELLQDVRV NQSRCEMKT VSPEAISPDL VKKIRDFHRK ILTLPEMMRM FSENLAHLE IDSGVITLDP QTASRSLVLS EDRKSVRYTR QKKSLPDSPL RFDGLPAVLG FPGFSSGRHR WQVDLQLGDG GGCTVGVAGE GVERRKGEMGL SAEDGVWAVI ISHQQCWAST SPGTDLPLSE IPRGVRVALD YEAGQVTLHN AQTQEPIFTF TASFSGKVFP FFAVWKKGSC LTLKG
<b>Source</b>	Mammalian cell
<b>Target Names</b>	TRIM15
<b>Protein Names</b>	Recommended name: Tripartite motif-containing protein 15 Alternative name(s): RING finger protein 93 Zinc finger protein 178 Zinc finger protein B7
<b>Expression Region</b>	1-465
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This protein is a member of the tripartite motif (TRIM) family. The TRIM motif includes three zinc-binding domains, a RING, a B-box type 1 and a B-box type 2, and a coiled-coil region. The protein localizes to the cytoplasm. Alternatively spliced transcript variants have been described, but their biological validity has not been determined.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.