



Recombinant Arabidopsis thaliana Protein ATAF2 (NAC081)

Product Code	CSB-EP887210DOA
Abbreviation	NAC081
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9C598
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MKSELNLPAG FRFHPTDEEL VKFYLCRKCA SEQISAPVIA EIDLKYNPW ELPEMSLYGE KEWYFFSPRD RKYPNGSRPN RAAGTGYWKA TGADKPIGKP KTLGIKKALV FYAGKAPKGI KTNWIMHEYR LANVDRSASV NKKNNLRLDD WVLCRIYNKK GTMEKYFPAD EKPRTTTMAE QSSSPFDTS STYPTLQEDD SSSSGGHGHV VSPDVLEVQS EPKWGELEDA LEAFDTSMFG SSMELLQPDA FVPQFLYQSD YFTSFQDPPE QKPFLNWSFA PQG
Source	E.coli
Target Names	NAC081
Protein Names	Recommended name: Protein ATAF2 Alternative name(s): NAC domain-containing protein 81 Short name= ANAC081
Expression Region	1-283
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.