



Recombinant Arabidopsis thaliana Probable WRKY transcription factor 9 (WRKY9)

Product Code	CSB-BP887247DOA
Abbreviation	WRKY9
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9C9F0
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MGFDFSTSKS KAKRQKRIEV RFASPLMGID LSLKLEAEEK KKEIEGSKHS RENKEDEEHD ASGDEDEQMV KEDEDDSSSL GLRTREEENE REELLQLLIQ MESVKEENTR LRKLVEQTL DYRHLEMKFP VIDKTKKMDL EMFLGVQGKR CVDITSKARK RGAERPSME REIGLSLSLE KKQKQEEESKE AVQSHHQRYN SSSLDMNMPR IISSSQGNRK ARVSVRARC TATMNDGCQW RKYGQKTAKG NPCPRAYYRC TVAPGCPVRK QVQRCLEDMS ILITTYEGTH NHPLPVGATA MASTASTSPF LLLDSSDNL HPSYYQTPQA IDSSLITYPQ NSSYNNRTIR SLNFDGPSRG DHVSSSQNRL NWMM
Source	Baculovirus
Target Names	WRKY9
Protein Names	Recommended name: Probable WRKY transcription factor 9 Alternative name(s): WRKY DNA-binding protein 9
Expression Region	1-374
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.