



Recombinant Arabidopsis thaliana E3 ubiquitin-protein ligase SINA-like 4 (At1g66650)

Product Code	CSB-EP887249DOA-B
Abbreviation	At1g66650
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9C9M0
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MTKLGRRNDG GGKSHRSSTK RQRRTSVSVD DPSPGEEEEK TLVVLTDSD SEEDDKPLGD VLRTCRKRRV SSPKSVTLPN SNVLECPNCF DPLKKPIFQC NNGHLACFLC CIKLLKRCFS CKLPIGDVRC RAMEKVIKAG LVSCSNAIYG CKQSTTYGNQ LQSHEKVCVF APCSCPIKDC NYIGFYKDLI NHFRATHKVS PGDINSFVFD RPFVIFGLDL DSSDKMVIFVE EKQGNLFVVQ GFIGSHGVYA TVSHIAPMVP EVRKFSCSLA RLRPYSTLRL GLEVKNIQKL RSQEEQPQED FLLIPSYMVN GDHMKMEISI GDKNDYTHI
Source	E.coli
Target Names	At1g66650
Protein Names	Recommended name: E3 ubiquitin-protein ligase SINA-like 4 EC= 6.3.2.- Alternative name(s): Seven in absentia-like protein 4
Expression Region	1-329
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.