



Recombinant Arabidopsis thaliana Putative clathrin assembly protein At1g68110 (At1g68110)

Product Code	CSB-BP887251DOA
Abbreviation	At1g68110
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9C9X5
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MKLWKRAAAA IKDRKSL LAV GFSRRNSSYR NADLEAAIIK ATSHDDSSVD YSNAHRVYKW IRSSPLNLKT LVYAISSRVN HTRSWIVALK SLMLLHG VLC CKVPSVVG EF RRLPFDLSDF SDGHSCLSKT WGFNVFVRTY FAFLLHHYSSF LSDQIHRLRG NNRRSLEKTS DSVIQELERI QKLQSLDMI LQIRPVADNM KKTLLILEAMD CLVIESINIY GRICGAVMKV LPLAGKSEAA TVLKIVNKTT SQGEDLIVYF EFCKGFGVSN AREIPQFVRI PEEEEVEAIEK MIDTVQE KPK LEKDEEKEDE KAMVVLEQPK KLQTIITDKW EIFEDDYRCF DRKDKWEIFE DEYHQNH LPL ITMNQPVYIT YTMPDLITF
Source	Baculovirus
Target Names	At1g68110
Protein Names	Recommended name: Putative clathrin assembly protein At1g68110
Expression Region	1-379
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.