



# Recombinant Mouse STAM-binding protein (Stampb)

<b>Product Code</b>	CSB-EP887415MO-B
<b>Abbreviation</b>	Stampb
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9CQ26
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSDHGDVSLP PQDRVRILSQ LGSVELNED IPPRRYYRSG VEIIRMASVY SEEGNIEHAF ILYNKYITLF IEKLPKHRDY KSIIPEKGD AVKKLKSVAF PKAEELKTEL LRRYTKEYEQ YKERKKKEEE ELARNIAIQQ ELEKEKQRVA QQKQKQLEQE QFHAFEEIMQ RQELEKERLK IVQEFQKQVDP GPCGPLLPLD EKPCVDVAPS SPFSPTQTPD CNTGMRPAKP PVVDRSLKPG ALSVIENVPT IEGLRHIVVP RNLCSFLQL ASANTAKGIE TCGVLCGKLM RNEFTITHVL IPRQNGGPDY CHTENEEEIF FMQDDLGLLT LGWIHTHTPTQ TAFLSSVDLH THCSYQMMLP ESIAIVCSPK FQETGFFKLT DYGLQEISTC RQKGFHPHGR DPPLFCDCSH VTVKDRIVTI TDLR
<b>Source</b>	E.coli
<b>Target Names</b>	Stampb
<b>Protein Names</b>	Recommended name: STAM-binding protein EC= 3.4.19.- Alternative name(s): Associated molecule with the SH3 domain of STAM
<b>Expression Region</b>	1-424
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	Cytokine-mediated signal transduction in the JAK-STAT cascade requires the involvement of adaptor molecules. One such signal-transducing adaptor molecule contains an SH3 domain that is required for induction of MYC and cell growth. This protein binds to the SH3 domain of the signal-transducing adaptor molecule, and plays a critical role in cytokine-mediated signaling for MYC induction and cell cycle progression. Multiple alternatively spliced transcript variants encoding the same protein isoform have been found for this gene.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.