



Recombinant Mouse Retinoic acid early-inducible protein 1-epsilon (Raet1e)

Product Code	CSB-MP887480MO
Abbreviation	Raet1e
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9CZQ6
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	LD DAHSLRCNLT IKDPTSADLP WCDVKCSVDE ITILHLNNIN KTMTSGDPGK MANATGKCLT QPLNDLCQEL RDKVSNTKVD THKTNGYPHL QVTMIYPQSQ GQTPSATWEF NISDSYFFTF YTENMSWRSA NDESGVIMNK WKDDGDLVQQ LKYFIPQCRQ KIDEFLKQSK EKPRSTSRSP SITQLTSTSP LPPPSHS
Source	Mammalian cell
Target Names	Raet1e
Protein Names	Recommended name: Retinoic acid early-inducible protein 1-epsilon Short name= RAE-1-epsilon
Expression Region	29-227
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.