



# Recombinant Mouse WAP four-disulfide core domain protein 2 (Wfdc2)

<b>Product Code</b>	CSB-EP887573MO-B
<b>Abbreviation</b>	Wfdc2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9DAU7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	TGTDA EKPGECPQLE PITDCVLECT LDKDCADNRK CCQAGCSSVC SKPNGPSEGE LSGTDTKLSE TGTTTQSAGL DHTTKPPGGQ VSTKPPAVTR EGLGVREKQG TCPSVDIPKL GLCEDQCQVD SQCSGNMKCC RNGCGKMACT TPKF
<b>Source</b>	E.coli
<b>Target Names</b>	Wfdc2
<b>Protein Names</b>	Recommended name: WAP four-disulfide core domain protein 2 Alternative name(s): WAP domain-containing protein HE4
<b>Expression Region</b>	26-174
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This gene encodes a protein that is a member of the WFDC domain family. The WFDC domain, or WAP Signature motif, contains eight cysteines forming four disulfide bonds at the core of the protein, and functions as a protease inhibitor in many family members. This gene is expressed in pulmonary epithelial cells, and was also found to be expressed in some ovarian cancers. The encoded protein is a small secretory protein, which may be involved in sperm maturation.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.