



Recombinant *Oryzias luzonensis* G2/mitotic-specific cyclin-B1 (ccnb1)

| | |
|--------------------------|--|
| Product Code | CSB-EP887615OCAG-B |
| Abbreviation | ccnb1 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9DG97 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Oryzias luzonensis</i> (Luzon ricefish) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MALRVTRNRL ASARAE LGGK TCSVAGPTQK PRAALGEIGN VAVINKDVTK KNVKTEVAKK TKVPAKAEKI EQPKAAVVPV KPAPVQVPA QPEPASPTPM ETSGCEPADL CQAFSDVILN TAIRDVDADD YNNPLLSEY VKDIYKYLRQ LEVEQSVKPN YLEGQEV TGN MRALLIDWLV QVSLKFRLQ ETMYMTVGII DRFLQDHPVP KKQLQLVGV T AMFLASKYEE MYPPEISDFA YVTDKAYTTA QIRDMEMTIL RVLKFQLGRP LPLQLRRAS KIYEVTAEQH TLAKYLLELS MVDYAMDHFP PSMVASAALA LTLKVL DAGE WDVTLQHYMA YTADTLTPVM AHIKNNVVKV NNGQTKHMTI KGKYFTSKQM RIATIFQLKS SVVKDLATQI S |
| Source | E.coli |
| Target Names | ccnb1 |
| Protein Names | Recommended name: G2/mitotic-specific cyclin-B1 |
| Expression Region | 1-401 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |