



# Recombinant Mouse WW domain-containing transcription regulator protein 1 (Wwtr1)

<b>Product Code</b>	CSB-EP887643MO-B
<b>Abbreviation</b>	Wwtr1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9EPK5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MNPSSVPHPL PPPGQQVIHV TQDLDTDLEA LFNSVMNPKP SSWRKKILPE SFFKEPDSGS HSRQSSTDSS GGHPGPRLAG GAQHVRSHSS PASLQLGTGA GAAGGPAQQH AHLRQQSYDV TDELPLPPGW EMTFTATGQR YFLNHIEKIT TWQDPRKVMN QPLNHVNLHP SITSTSVQR SMAVSQPNLA MNHQHQVVA TSLSPQNHPT QNQPGLMSV PNALTTQQQQ QQKLRLQRIQ MERERIRMRQ EELMRQEAAL CRQLPMETET MAPVNTPAMS TDMRSVTNSS SDPFLNGGPY HSREQSTDSG LGLGCYSVPT TPEDFLSNMD EMDTGENSGQ TPMTVNPQQT RFPDFLDCLP GTNVDLGTLE SEDLIPLFND VESALNKSEP FLTWL
<b>Source</b>	E.coli
<b>Target Names</b>	Wwtr1
<b>Protein Names</b>	Recommended name: WW domain-containing transcription regulator protein 1 Alternative name(s): Transcriptional coactivator with PDZ-binding motif
<b>Expression Region</b>	1-395
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.