



# Recombinant Dickeya dadantii Phosphatidylserine decarboxylase proenzyme (psd)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-YP887682DKD   |
| <b>Abbreviation</b>      | psd   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q9EV04  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Dickeya dadantii (strain 3937) (Erwinia chrysanthemi (strain 3937))   |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MLDRIKIALQ HLLPKVWLTQ LAGWGADRQA GMLTKLVIDL FARIYKVNMQ<br>EAQQPDTASY RSFNDFVVRP LKPGIRPVDL LPNRLVFPAD GAISQLGAID<br>DDRILQAKQH DYTLEALLAG NYIISDLFRD GLFVTTYLSP RDYHRVHMPC<br>DGILRDMIYV PGDLFSVNPL TAANVPNLFA RNERVICLFD TPFGPMVQIL<br>VGATIVGSIE TVWAGVTPP REGIIKRWAY PMELEGAVIL EKGDEMGRFK<br>LG   |
| <b>Source</b>            | Yeast   |
| <b>Target Names</b>      | psd   |
| <b>Protein Names</b>     | Recommended name: Phosphatidylserine decarboxylase proenzyme EC=4.1.1.65 Cleaved into the following 2 chains: 1. Phosphatidylserine decarboxylase alpha chain 2. Phosphatidylserine decarboxylase beta chain  |
| <b>Expression Region</b> | 1-252   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |