



Recombinant Arabidopsis thaliana B3 domain-containing protein At5g38500 (At5g38500)

Product Code	CSB-EP887739DOA-B
Abbreviation	At5g38500
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FFX1
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MDSGKDMWSR LCLLAETVVM AAEEEEQRRR LLAEKREDSK SQKKTVSEED DSEKRFLSHV PRKKRSSLVK RQQKPNGVST SSSSLPDLNQ IPIDYETETK QNPSFIERLV CDEEQRVKKG KSRIWEEEE EADEDESEKRL FEKNLMKFVG HSQQQKQFET LNGASSSSSF LNLRCYEASL FLDYNTVESE KTETKVL PNP NYQSSSPSSC LTENDTSRKR RAVEQRKSGK VKKVKVSPLP RLSTETPEWV FQAMGHMNAD AETPKLIFER TLFKSDVNSN LSRLIPFQK LIRNDFLTPE ECRAMQEDKD KDDDISVGT ILVCQAKQED EDKDDDEDIGA GTILVNQRFK MWGLRFKIWG MEKDSGHGTL NYILNWDWND VVKGNLSKAG DNIGLWTFRC RGVLCFALDT W
Source	E.coli
Target Names	At5g38500
Protein Names	Recommended name: B3 domain-containing protein At5g38500
Expression Region	1-411
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.