



# Recombinant Arabidopsis thaliana B3 domain-containing protein At5g38500 (At5g38500)

<b>Product Code</b>	CSB-MP887739DOA
<b>Abbreviation</b>	At5g38500
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9FFX1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDSGKDMWSR LCLLAETVVM AAEEEEQRRR LLAEKREDSK SQKKTVSEED DSEKRFLSHV PRKKRSSLVK RQQKPNGVST SSSSLPDLNQ IPIDYETETK QNPSFIERLV CDEEQRVKKG KSRIWEEEE EADEDSEKRL FEKNLMKFVG HSQQQQKFET LNGASSSSSF LNLRCYEASL FLDYNTVESE KTETKVL PNP NYQSSSPSSC LTENDTSRKR RAVEQRKSGK VKKVKVSPLP RLSTETPEWV FQAMGHMNAD AETPKLIFER TLFKSDVNSN LSRLIPFQK LIRNDFLTPE ECRAMQEDKD KDDDISVGT ILVCQAKQED EDKDDDEDIGA GTILVNQRFK MWGLRFKIWG MEKDSGHGTL NYILNWDWND VVKGNLSKAG DNIGLWTFRC RGVLCFALDT W
<b>Source</b>	Mammalian cell
<b>Target Names</b>	At5g38500
<b>Protein Names</b>	Recommended name: B3 domain-containing protein At5g38500
<b>Expression Region</b>	1-411
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.