



Recombinant Arabidopsis thaliana Putative F-box protein At5g50220 (At5g50220)

Product Code	CSB-MP887748DOA
Abbreviation	At5g50220
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FGS3
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MKRGSEEIKE ETSPLPAKHF QRREVSIAED IGIPIDLMVE ILKKLPAKSL IKFQCVSKQW SSIIGSSRDF IDSIVTRSL S QPSRDILISF STLTNSLKQ ISSSFPLRTL DILTKNQSYT EAAIYNPTTR QLSLSPETTA GHSHVSTSFL GYDPFKNQYK VICLDNYKRR CCHVFTLGDA IRKWRKIQYN FGLYFPLLPP VCIKGTIYYQ AKQYGSTYVL LCFDWISEKF DQVEAPKTMM DHRYTLINYQ GKLGFMCCQN RVEIWVMKND EKKQEWKIF FYEMAGFEKW HIARATPSGE IVFVNRLLLS CQTLVYYYYG PKRNSMRRVE VEGTKYRRKH LVHICPVPDH VENTMRL
Source	Mammalian cell
Target Names	At5g50220
Protein Names	Recommended name: Putative F-box protein At5g50220
Expression Region	1-357
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.