



Recombinant Arabidopsis thaliana Probable WRKY transcription factor 49 (WRKY49)

Product Code	CSB-EP887759DOA-B
Abbreviation	WRKY49
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FHR7
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MEEEGYQWAR RCGNNAVEDP FVYEPPLFFL PQDQHMHMGL MPNEDFIANK FVTSTLYSGP RIQDIANALA LVEPLTHPVR EISKSTVPLL ERSTLSKVDR YTLKVKNNNSN GMCDDGYKWR KYGQKSIKNS PNPRSYYKCT NPICNAKKQV ERSIDESNTY IITYEGFHFH YTYPPFLPDK TRQWPNKKTK IHKHNAQDMN KKSQTQEEK EAQLGELTNQ NHPVNKAQEN TPANLEEGLF FPVDQCRPQQ GLLLEDVVAPA MKNIPTRDSV LTAS
Source	E.coli
Target Names	WRKY49
Protein Names	Recommended name: Probable WRKY transcription factor 49 Alternative name(s): WRKY DNA-binding protein 49
Expression Region	1-274
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.