



# Recombinant Arabidopsis thaliana GDSL esterase/lipase At5g45910 (At5g45910)

<b>Product Code</b>	CSB-EP887771DOA-B
<b>Abbreviation</b>	At5g45910
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9FJ45
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	L PMRPTLKYES IFNFGDSLSD TGNFLLSGDV DSPNIGRLPY GQTFNFRSTG RCSDGRLIID FIAEASGLPY IPPYLQSLRT NDSVDFKRGGA NFAVAGATAN EFSFFKNRGL SVTLLTNKTL DIQLDWFKKL KPSLCKTKPE CEQYFRKSLF LVGEIGGNDY NYPLLAFRSF KHAMDLVPFV INKIMDVTSALIEEGAMTLI VPGNLPIGCS AALLERFNDN SGWLYDSRNQ CYMPLNNAK LHNKLLKGL AALRKKYPYA KIIYADYYSS AMQFFNSPSK YGFTGSVLKA CCGGGDGRYN VQPNVRCGEK GSTTCEDPST YANWDGIHLT EAAYRHATG LISGRFTMPT YN
<b>Source</b>	E.coli
<b>Target Names</b>	At5g45910
<b>Protein Names</b>	Recommended name: GDSL esterase/lipase At5g45910 EC= 3.1.1.- Alternative name(s): Extracellular lipase At5g45910
<b>Expression Region</b>	20-372
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.