



Recombinant Arabidopsis thaliana Pentatricopeptide repeat-containing protein At5g48730, chloroplastic (At5g48730)

Product Code	CSB-BP887785DOA
Abbreviation	At5g48730
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FKC3
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	CI SISPREPNYA ITSDKSNNTS LSLRETRQSK WLINAEDVNE RDSKEIKEDK NTKIASRKAI SIILRREATK SIIKKKGGSK KLLPRTVLES LHERITALRW ESAIQVFELL REQLWYKPNV GIYVKLIVML GKCKQPEKAH ELFQEMINEG CVVNHEVYTA LVSAYSRSGR FDAAFLLER MKSSHNCQPD VHTYSILIKS FLQVFAFDKV QDLLSDMRRQ GIRPNTITYN TLIDAYGKAK MFVEMESTLI QMLGEDDCKP DSWTMNSTLR AFGGNGQIEM MENCYEKFQS SGIEPNIRTF NILDSYGKS GNYKKMSAVM EYMQKYHYSW TIVTYNVVID AFGRAGDLKQ MEYLFRLMQS ERIFPSCVTL CSLVRAYGRA SKADKIGGVL RFIENSDIRL DLVFFNCLVD AYGRMEKFAE MKGVLELMEK KGFKPKDKITY RTMVKAYRIS GMTTHVKELH GVVESVGAEQ VVVKKPDF
Source	Baculovirus
Target Names	At5g48730
Protein Names	Recommended name: Pentatricopeptide repeat-containing protein At5g48730, chloroplastic
Expression Region	29-508
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.