



Recombinant Arabidopsis thaliana Putative E3 ubiquitin-protein ligase SINA-like 9 (At5g37910)

| | |
|--------------------------|---|
| Product Code | CSB-BP887786DOA |
| Abbreviation | At5g37910 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9FKD5 |
| Product Type | Recombinant Protein |
| Immunogen Species | Arabidopsis thaliana (Mouse-ear cress) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | MVLASISEAL ISQGDGGERV AKRQRSAIVL LDLDILDCPI CCEALTSPIF QCDNGHLACG SCCPKLSNKC PACTLPVGH SRSRAMESVLE SILIPCPNVR FGCTKSFFYG KESAHEKECI FSQCSCPSSV CDYTGSYKDL YAHYKLTHST NIFWNIKRFR CANFFTT SML ISDKILIKRV HEKKLLAVQ CFREPCGVYV TVSFIAPSAP EVGEFSYQLS YNVDGHTVTY ESPEVKRVCK VSIETPQENF MLIPHSLLRG DLLEMQVFII ENV DQE |
| Source | Baculovirus |
| Target Names | At5g37910 |
| Protein Names | Recommended name: Putative E3 ubiquitin-protein ligase SINA-like 9 EC=6.3.2.- Alternative name(s): Seven in absentia-like protein 9 |
| Expression Region | 1-276 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |