



Recombinant Arabidopsis thaliana Putative E3 ubiquitin-protein ligase SINA-like 9 (At5g37910)

Product Code	CSB-MP887786DOA
Abbreviation	At5g37910
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FKD5
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MVLASISEAL ISQGDGGERV AKRQRSAIVL LDLDILDCPI CCEALTSPIF QCDNGHLACG SCCPKLSNKC PACTLPVGH SRSRAMESVLE SILIPCPNVR FGCTKSFFYG KESAHEKECI FSQCSCPSSV CDYTGSYKDL YAHYKLTHST NIFWNIKRFR CANFFTT SML ISDKILIKRV HEKKLLAVQ CFREPCGVYV TVSFIAPSAP EVGEFSYQLS YNVDGHTVTY ESPEVKRVCK VSIETPQENF MLIPHSLLRG DLLEMQVFII ENV DQE
Source	Mammalian cell
Target Names	At5g37910
Protein Names	Recommended name: Putative E3 ubiquitin-protein ligase SINA-like 9 EC= 6.3.2.- Alternative name(s): Seven in absentia-like protein 9
Expression Region	1-276
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.