



# Recombinant Arabidopsis thaliana Probable aminotransferase TAT2 (At5g53970)

<b>Product Code</b>	CSB-MP887820DOA
<b>Abbreviation</b>	At5g53970
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9FN30
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MENGATTTST ITIKGILSLL MESITTEEDE GGKRVISLGM GDPTLYSCFR TTQVSLQAVS DSLLSNKFHG YSPTVGLPQA RRAIAEYLSR DLPYKLSQDD VFITSGCTQA IDVALSMLAR PRANILLPRP GFPIYELCAK FRHLEVRYVD LLPENGWEID LDAVEALADE NTVALVVINP GNPCGNVYSY QHLMKIAESA KKGFLVIAD EVYGHAFGS KPFVPMGVFG SIVPVLTLGS LSKRWIVPGW RLGWFVTTDP SGSFKDPKII ERFKDYFDIL GGPATFIQAA VPTILEQTDE SFFKKTLSL KNSSDICCDW IKEIPCIDSS HRPEGSMAMM VKLNLSLLED VSDDIDFCFK LAREESVILL PGTAVGLKNW LRITFAADAT SIEEAFKRIK CFYLRHAKTQ YPTI
<b>Source</b>	Mammalian cell
<b>Target Names</b>	At5g53970
<b>Protein Names</b>	Recommended name: Probable aminotransferase TAT2 EC= 2.6.1.- Alternative name(s): Tyrosine aminotransferase 2
<b>Expression Region</b>	1-414
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.