



Recombinant Arabidopsis thaliana BTB/POZ domain-containing protein At5g41330 (At5g41330)

Product Code	CSB-BP887821DOA
Abbreviation	At5g41330
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FN67
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MNFPTIPHES NVVSINVGGR IFQTTKQTLS LAGTDSLLSQ LATETTRFVD RDPDLFSVLL YILRTGNLPA RSRAFDVRL IDESRYYGIE SFLIDSLSNS SQFEPDLRR SRILQLNGRD SPSSISPTVI GGGLHVAHGS KITSFDWSLR QKSTVLTNFS AVDSLLEISP GVLAAGATDF PGLQIIDLDN GGFVRTTLNW ENVTRSSSTV QAIGSSHEFL FTSFESSRRN SNSIMVYDLS SLLPVSEIDH CEIYGADIDS AIPSTKLKWI QSCNLLMVSG SHTSPSGVNG HIRFWDVRSR NMVWEIKEAQ DCFSDVTVSD NLSAVFKVGV TSGEVFYADL RSLGTKDPWV CLGEEKRNL NERRGVGSKI ESYGNHVFCS SKGSGIELWS EVITGLVGNA SRDVLEERVF RKNSLGKLAD SGENKITGLA FGGNRMFVTR KDQQSVQVWQ SPSRGISI
Source	Baculovirus
Target Names	At5g41330
Protein Names	Recommended name: BTB/POZ domain-containing protein At5g41330
Expression Region	1-458
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.