



Recombinant Arabidopsis thaliana F-box/FBD/LRR-repeat protein At5g22660 (At5g22660)

Product Code	CSB-MP887825DOA
Abbreviation	At5g22660
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FNJ5
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MSNQGAIARRS GEDRISSLPD HLLSQILSNL PTENAVTTSI LSTRWKDLWL STPVLDIDID AFDDATTFIS FATRFLDSFK DSCLHKLQIS FQMEAVDMWT IIPWIEDAVK RRIQHLEVDS RIDHMIDTLP LTVYLSESLV SLRLHLVMLH RFV FVSLPNL KVMHLEENIY SYAETMEKFI SSCPVLEDLT VVRNVDEATE KVL RVSSQSL NSLKLVIDSS KCWYNDDSDD WKVVIDAPQL VYLSLKDDQS VSFVINNLCS SAKADIKVSF NVSDIWDLEE SFERSNVGKF LTGLSSLRDM TISGTTLKII CHYLKHEPMP QFRNMTRLHA KFYVCDLEML PCVLESCP NL KSLVLKLGKE MENEEISLSS SVPKCLQSSL ENVEIRPNY GSGEEMKLSK YFLENSLVLK KFKLCRDCHS EEQESLVVRE LMTFQRCSSA CEINVVRFQR
Source	Mammalian cell
Target Names	At5g22660
Protein Names	Recommended name: F-box/FBD/LRR-repeat protein At5g22660
Expression Region	1-450
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.