



Recombinant Arabidopsis thaliana Putative F-box protein At1g47765 (At1g47765)

Product Code	CSB-BP887874DOA
Abbreviation	At1g47765
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FZF3
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MEQQKQKKRK VVSKSKRTQS KSASSLPLDL TSEILLRLPE KSIARFRCVS KLWLSITDP YFINLFETRS PRPSLLVCFI ENDKLFVSSI PQHLHSLQNS KRSYSSSQLF IVIWNQYQKD VGLICFKTSK MPIVWNPSKR QLITLPIPRL SWNNIIVFLG YDPVEGKHKV MCLPFRRSSD VCQVLTGPA QEFSWITVKT YHKHCSDYQS SGRCIKGVVY YIAQVYHTHA WVL MCFD VRS EKFDMIKLHA DIYREILITY EGRIACVEKR TTKDDYIALC ILEDAKKDKW SSKDILAPFG HFDERLRTFF QLK GCTHDGE FIFVSSTFRK MDYILFFDPV KKTFRRFELK EIADDQARVN NGDPYGPIFA FCAFLDHVES QMSLQ
Source	Baculovirus
Target Names	At1g47765
Protein Names	Recommended name: Putative F-box protein At1g47765
Expression Region	1-385
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.