



Recombinant Anemonia sulcata GFP-like non-fluorescent chromoprotein FP595

Product Code	CSB-BP887932AKE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9GZ28
Product Type	Recombinant Protein
Immunogen Species	Anemonia sulcata (Mediterranean snakelocks sea anemone)
Purity	>85% (SDS-PAGE)
Sequence	MASFLKKTMP FKTTIEGTVN GHYFKCTGKG EGNPFEGTQE MKIEVIEGGP LPFAFHILST SC
Source	Baculovirus
Protein Names	Recommended name: GFP-like non-fluorescent chromoprotein FP595 Alternative name(s): asFP595 Cleaved into the following 2 chains: 1. GFP-like non-fluorescent chromoprotein FP595 chain 1 2. GFP-like non-fluorescent chromoprotein FP595 chain
Expression Region	1-62
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.