



Recombinant Human Cleavage stimulation factor subunit 2 tau variant (CSTF2T)

Product Code	CSB-EP887946HU-B
Abbreviation	CSTF2T
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H0L4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSSLAVRDPA MDRSLRSV FV GNIPYEATEE QLKDIFSEVG SVVSFRLVYD RETGKPKGYG FCEYQDQETA LSAMRNLNGR EFSGRALRVD NAASEKNKEE LKSLGPAAPI IDSPYGDPID PEDAPESITR AVASLPPEQM FELMKQMKLC VQNSHQEARN MLLQNPQLAY ALLQAQVVMR IMDPEIALKI LHRKIHVTPL IPGKSQSVSV SGP GPGPGPG LCPGPNVLLN QQNPPAPQPQ HLARRPVKDI PPLMQTPIQG GIPAPGPIPA AVPGAGPGSL TPGGAMQPQL GMPGVGPVPL ERGQVQMSDP RAPIRGPVT PGGLPPRGLL GDAPNDPRGG TLLSVTGEVE PRGYLGPPHQ GPPMHASGH DTRGPSSHEM RGGPLGDPRL LIGEP RGP MI DQRGLPMDGR GGRDSRAMET RAMETE VLET RVMERRGMET CAMETRGMEA RGMDARGLEM RGPVPSSRGP MTGGIQGPGP INIGAGGPPQ GPRQVPGISG VGNPGAGMQG TGIQGTGMQG AGIQGGGMQG AGIQGVSIQG GGIQGGGIQG ASKQGGSQPS SFSPGQSQVT PQDQEKAALI MQVLQLTADQ IAMLPP EQRQ SILILKEIQ KSTGAS
Source	E.coli
Target Names	CSTF2T
Protein Names	Recommended name: Cleavage stimulation factor subunit 2 tau variant Alternative name(s): CF-1 64 kDa subunit tau variant Cleavage stimulation factor 64 kDa subunit tau variant Short name= CSTF 64 kDa subunit tau variant TauCstF-64
Expression Region	1-616
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.